



One Step Test for NGAL

(Colloidal Gold)

User Manual

Cat.# CG1010

INTENDED USE

One Step Test for NGAL (Colloidal Gold) is intended for *in vitro* quantitative determination of neutrophils gelatinase associated lipocalin (NGAL) in serum and urine. This test is used as an aid in the early diagnosis of acute kidney injury (AKI), risk classification and treatment monitoring.

SUMMARY

The inclining incidence of chronic kidney disease which has led to high mortality and immense medical burden over the past decades has become a distressing concern in epidemiology. Unfortunately, the number of biomarkers that allow the monitoring of chronic kidney disease (CKD) is limited. NGAL is an emerging biomarker which has been shown to be able to aid the diagnosis of kidney injuries.

The evidence for the role of NGAL measurements in a variety of clinical situations leading to AKI (cardiac surgery, kidney transplantation, contrast nephropathy, haemolytic uraemic syndrome and in the intensive care setting) or to CKD (lupus nephritis, glomerulonephritides, obstruction, dysplasia, polycystic kidney disease, IgA nephropathy) is explored. The emerging utility of standardized clinical platforms for reliable measurement of NGAL in plasma (Triage NGAL Device; Biosite Incorporated) and urine (ARCHITECT analyzer; Abbott Diagnostics) is also discussed. It will be important in future studies to validate the sensitivity and specificity of NGAL concentration measurements in clinical samples from large cohorts and from multiple clinical situations. Such studies will be facilitated by the anticipated widespread availability of standardized commercial tools in the near future.

PRINCIPLE

The test uses an anti-human NGAL polyclonal antibody conjugated with colloidal gold and an anti-human NGAL monoclonal antibody coated on the test line. After the sample has been applied to the test strip, the gold-labelled anti-human NGAL polyclonal antibody binds with the NGAL in sample and forms a marked antigen-antibody complex. This complex moves to the test card detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human NGAL monoclonal antibody resulting in a purplish red streak appears on the test line. The color intensity of the test line increases in proportion to the amount of NGAL in sample.

Then insert test card into FIA8000 Quantitative Immunoassay Analyzer (hereinafter referred to as FIA8000), the concentration of NGAL in sample will be measured and displayed on the screen. The value will be stored in FIA8000 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

A kit contains:

- 1. Getein NGAL test card in a sealed pouch with desiccant 25
- 2. Disposable pipet 25
- 3. User manual 1
- 4. SD card 1
- 5. Sample diluent 25

A test card consists of:

A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (a colloidal gold-labelled anti-human NGAL polyclonal antibody is coated at the border of the nitrocellulose membrane and sample pad, the test line is coated with an anti-human NGAL monoclonal antibody, and the control line is coated with rabbit anti-goat IgG antibody), absorbent paper and liner.

Sample diluent:

Phosphate buffered saline, proteins, detergent, preservative, stabilizer.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

FIA8000 Quantitative Immunoassay Analyzer

STORAGE AND STABILITY

Store the test card at 4~30°C with a valid period of 24 months. Use the test card within 1 hour once the foil pouch is opened.

PRECAUTIONS

- 1. For *in vitro* diagnostic use only.
- 2. For professional use only.
- 3. Do not use the kit beyond the expiration date.
- 4. Do not use the test card if the foil pouch is damaged.
- 5. Do not open pouches until ready to perform the test.
- 6. Do not reuse the test card.
- 7. Do not reuse the pipet.
- 8. Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
- 9. Carefully read and follow the manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

- 1. This test can be used for *serum and urine samples*. Samples should be free of hemolysis.
- 2. If testing will be delayed, serum sample may be stored up to 7 days at 2~8°C or stored at -20°C for 6 months before testing .
- 3. Urine sample can be preserved at room temperature for 4 hours, please test it as soon as possible. If testing will be delayed, urine sample may be stored up to 3 days at 2~8°C before testing.
- 4. Samples should be brought to room temperature before testing.
- 5. Do not use frozen urine samples.
- 6. Do not use heat-inactivated samples.
- 7. SAMPLE VOLUME: **10 µl**.

TEST PROCEDURE

- 1. Collect specimen according to user manual.
- 2. Test card, sample should be brought to room temperature

before testing.

3. Confirm SD card lot No. in accordance with test kit lot No.. Perform "QC (SD)" calibration when necessary (Details refer to 8.2.1 of FIA8000 User Manual).
4. On the main interface of FIA8000, press "ENT" button to enter testing interface.
5. Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
6. Put the test card on a clean table, horizontally placed.
7. Using sample transfer pipette, deliver **10 µl** of sample into one tube of sample diluent, mix gently and thoroughly. Then drop **100 µl** of sample mixture (or 4 drops of sample mixture when using disposable pipet) into the sample port on the test card.
8. **Reaction time: 3 minutes.** Insert the test card into FIA8000 and press "ENT" button after reaction time is elapsed. The result will be shown on the screen and printed automatically.

Notes:

1. It is required to perform "QC (SD)" calibration when using a new batch of kits.
2. It is suggested to calibrate once for one batch of kits.
3. Make sure the test card insertion is correct and complete.

TEST RESULTS

Valid: When a purplish-red band appears at the control area (C), use FIA8000 to analyze the test card and get the result.

Invalid: If no colored band appears in the control area (C), the test result is invalid. The test should be repeated and if the same situation happened again, please stop using this batch of products and contact your supplier.

EXPECTED VALUE

The expected normal value for NGAL was determined by testing samples from 319 apparently healthy individuals. The 95th percentile of the concentration for NGAL in serum is 200 ng/ml. The 95th percentile of the concentration for NGAL in urine is 100 ng/ml. (The probability that value of a normal person with serum below 200 ng/ml is 95%, the probability that value of a normal person with urine below 100 ng/ml is 95%.) It is recommended that each laboratory establish its own

expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range	50–5000 ng/ml
Lower Detection Limit	≤50 ng/ml
Within-Run Precision	≤10%
Between-Run Precision	≤15%
Method Comparison:	

The assay was compared with HITACHI 7600 analyzer and Beijing Strong NGAL Test Kits with 200 urine and serum samples. The correlation coefficient (r) is 0.980, 0.989 respectively.

LIMITATIONS

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
2. Samples containing interferences may influence the results. The table below listed the maximum allowance of these potential interferences.

Interferent	Creatinine	Glucose	Urea
Concentration (Max)	10 g/L	10 g/L	100 g/L
Interferent	Hemoglobin	Triglyceride	Bilirubin
Concentration (Max)	5 g/L	25 g/L	0.1 g/L













REFERENCES

1. Wasilewska A, Taranta-Janusz K, Dębek W, et al. KIM-1 and NGAL: new markers of obstructive nephropathy. *Pediatr Nephrol.* 2011, 26(4):579-586.
2. Clerico A, Galli C, Fortunato A, et al. Neutrophil gelatinase-associated lipocalin (NGAL) as biomarker of acute kidney injury: a review of the laboratory characteristics and clinical evidences. *Clin Chem Lab Med.* 2012, 50(9):1505-1517.
3. Shemin D, Dworkin LD. Neutrophil gelatinase-associated lipocalin (NGAL) as a biomarker for early acute kidney injury. *Crit Care Clin.* 2011, 27(2):379-389.
4. EN ISO 18113-1:2009 *In vitro* diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements.
5. EN ISO 18113-2:2009 *In vitro* diagnostic medical devices -

Information supplied by the manufacturer (labelling) - Part 2: *In vitro* diagnostic reagents for professional use (ISO 18113-2:2009).

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on One Step Test for NGAL (Colloidal Gold) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2007.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		<i>In vitro</i> diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing One Step Test for NGAL (Colloidal Gold). Please read this user manual carefully before operating to ensure proper use.

Version: WCG14-DL-S-03



Getein Biotech, Inc.

Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China

Tel: +86-25-68568508

Fax: +86-25-68568500

E-mail: tech@getein.com.cn

overseas@getein.com.cn

Website: www.bio-GP.com.cn